

## Prody® Abdominal Binder II with Drain Fasteners

The Abdominal Binder II with drain fasteners is the newest version of the original Abdominal Binder, with improved easy grip fasteners making drain care easier. It securely supports and organizes drainage bulbs and allows patients to move more comfortably while promoting drain function.

### Prody® Abdominal Binder II with Drain Fasteners

| Panels | Width | Stretch   | Unit | Model #    |
|--------|-------|-----------|------|------------|
| 2      | 6"    | Up to 60" | Each | PPABD-601  |
| 3      | 9"    | Up to 45" | Each | PPABD-901  |
| 3      | 9"    | 46" – 62" | Each | PPABD-902  |
| 3      | 9"    | 55" – 72" | Each | PPABD-903  |
| 4      | 12"   | Up to 45" | Each | PPABD-1201 |
| 4      | 12"   | 46" – 62" | Each | PPABD-1202 |
| 4      | 12"   | 55" – 72" | Each | PPABD-1203 |
| 4      | 12"   | 72" – 96" | Each | PPABD-1204 |



US Patent #11,090,184

## Benefits and Advantages

- New improved easy grip fasteners that are easier to grip for drain care. Drainage Bulbs are visible to determine need to empty, yet not visible to others under clothing
- Secures up to four drainage bulbs eliminating the need for safety pins
- Maximizes drain function, decreases tugging, tangling or potential dislodgement
- Paneled elastic construction with plush backing
- Stitching between horizontal panels allows binder to conform to body shape
- Decreases complications and improves outcomes

