



Promedics Veterinary Products

Prody® II Post-Surgical Drain and Catheter Care Products

Exclusively Manufactured by:

NYOrtho

800-699-2358 www.nyorthousa.com



US Patent # 11,090,184



Prody® II Veterinary 2" Drain and Catheter Holder

Size:	Unit:	Model#	Case:
Universal	Each	8251	50

Prody® II Veterinary 6" Drain and Catheter Holder

Panels:	Height:	Fits:	Unit:	Model#
2	6"	30" - 45"	Each	8250-3045
2	6"	45" - 60"	Each	8250-4560

FEATURES

Drain and Catheter Security

- Designed with individual fasteners to secure post-surgical drains, catheters, and enteral feeding tubes.
- New easy-grip fasteners make drain and catheter care easier.
- Prevents tangling or accidental dislodgement of drains, catheters, or tubes.
- Secures up to four drainage bulbs.

Comfort, Fit and Outcomes

- Allows patients to move freely and comfortably.
- Wraps securely around the chest or abdomen.
- Available in a 2-inch universal width that stretches up to 54 inches and can be trimmed, or in a 6-inch width with two sizes for animals over 50 pounds.
- Secure line management helps decrease complications and improve outcomes.

CARE INSTRUCTIONS

- Before washing, securely attach the blue and white hook tabs onto the elastic fabric. Hand wash in cold water.
- Hang and allow to air dry completely.

For more information visit our website www.promedicsveterinaryproducts.com

To order: contact your distribution partner, email prody@nyorthousa.com or call **800.699.2358**



US Patent # 11,090,184